

### UNIVERSITI TUN HUSSEIN ONN MALAYSIA

# FINAL EXAMINATION SEMESTER I SESSION 2015/2016

**COURSE NAME** 

MEDICAL MEASUREMENT

**SYSTEM** 

COURSE CODE

: BEU41003

**PROGRAMME** 

BACHELOR OF ELECTRONIC

ENGINEERING WITH HONOURS

**EXAMINATION DATE** 

DECEMBER 2015 / JANUARY 2016

**DURATION** 

: 3 HOURS

:

INSTRUCTION

1. WRITE ALL THE

**ANSWERS** IN THIS

**QUESTIONS BOOKLET** 

2. ANSWER ALL QUESTIONS

THIS QUESTION PAPER CONSISTS OF THIRTEEN (13) PAGES

- Q1 (a) The instrumentation amplifier is a key design component to almost all biopotential measurement.
  - (i) Discuss the reasonable usage of the amplifier in almost all medical measurement system. (4 marks)

(ii) Derive the equations for the gain of Op-amp A<sub>3</sub> instrumentation amplifier in **Figure Q1(a)(ii)** and calculate the overall circuit gain. Given that:

$$R_2 = 50 \text{k}\Omega$$
,  $R_1 = 0.5 \text{k}\Omega$  and  $R_3 = R_4 = 50 \text{k}\Omega$ .

The output voltage from Op-amp  $A_1$  is  $V_3$  and the output voltage from Op-amp  $A_2$  is  $V_4$ .

(10 marks)




FIGURE Q1(a)(ii)

(b)	(i)	Explain the usage of isolation amplifiers in medical measurement system.  (3 marks)
	(ii)	Draw the block diagram of the basic design of an isolation amplifier and explain the function of each block.  (8 marks)

BEU 41003

Q2	(a)	Explain the principle of photoplethysmogram (PPG) and illustrate the PPG signal						
			(7 marks)					

(b) Draw and explain the **TWO (2)** modes of photoplethysmogram (PPG).

(8 marks)

BEU 41003


	Design the block diagram of photoplethysmogram (PPG).	(c)
(5 marks)		

(d) Explain briefly Pulse Transit Time (PTT).

(5 marks)

**Q3** (a) **Figure Q3(a)** shows the oscilloscope display of the blood pressure measurement using oscillometric sensor. From the result shown, calculate the diastolic pressure (mmHg) and systolic pressure (mmHg). Given that the VOLT/DIV for both channels (CH1 and CH2) are set at 2 V/div while the TIME/DIV are set at 5 Sec/DIV.

(6 marks)






FIGURE Q3(a)

BEU 41003

#### BEU 41003


(b)	Using the diastolic and	l systolic	pressure	measured	in Fig	gure Q3(a),	calculate	the
	mean arterial pressure.							
							(2 mai	rks)

- (c) We would like to monitor the heart rate of a subject while exercising. Exercise causes three main difficulties for measurement of the Electrocardiogram (ECG):
  - The movement of artifacts.
  - Sweat on the electrodes.
  - Artifacts from the Electromyogram (EMG).
  - (i) Explain the definition of artifacts.

(2 marks)

(ii) Determine the type of artifact that affect the ECG waveform in Figure Q3(c)(ii) and describe the possible causes. Explain if a filter can be used to rectify the signal?


(3 marks)



#### FIGURE Q3(c)(ii)

(iii) Set-up the connection of LEAD I, LEAD II and LEAD III for the ECG measurement in Figure Q3(c)(iii).

(6 marks)



#### FIGURE Q3(c)(iii)

(iv) Label the P, Q, R, S, T events in **Figure Q3(c)(iv)**. Estimate the heart rate (beats/min) based on the corresponding diagram.

(6 marks)



#### FIGURE Q3(c)(iv)

Q4 (a) (i) Explain the principles of respiratory ventilation detection. Use diagram if necessary.

(5 marks)

#### BEU 41003

(ii) Illustrate and compare the difference between measured waveforms of respiratory while relax with the respiratory while exercise.

(7 marks)

(b) (i) In the process of your heart pumps blood through your body, you can feel a pulsing in some of the blood vessels close to the skin's surface, such as in your wrist, neck, or upper arm. Demonstrate the counting of your pulse rate when your heart is beating.

(5 marks)

BEU 41003

(ii) Propose a complete method for transferring circuit principles of pulse rate measurement for heart rate analysis and measurement purposes.

(8 marks)

**END OF QUESTION -**